




Silverline®

Antimicrobial Peritoneal Catheters



Antimicrobial EVD Catheters

Silverline Catheters are impregnated with 2% of nanoparticles of silver and an insoluble silver salt. Silver ions are a strong antiseptic. Because of the nanometer size of the particles the specific surface is large, which results in the release of ions in antimicrobial concentration on the surface of the catheter.


The silver salt starts to release silver ions immediately when it gets in contact with body fluid. The release of silver ions from the silver salt stops after some hours. At that time the silver particles begin to release silver ions.

The combination of a fast acting substance and a slow acting substance results in an optimal short term and long term effect.

In three different shapes

Silverline Antimicrobial Ventricular Catheter 8F

Order Number EVD 30.010.02
Material Polyurethane
Outside Diameter 2.7 mm
Inside Diameter 1.6 mm
Length 270 mm



Silverline Antimicrobial Ventricular Catheter with Hollow Stylet 8F 'Heese'

Order Number EVD 30.012.02
Material Polyurethane
Outside Diameter 2.7 mm
Inside Diameter 1.6 mm
Length 200 mm



Silverline Antimicrobial Short Ventricular Catheter 8F

Order Number EVD 30.011.02
Material Polyurethane
Outside Diameter 2.7 mm
Inside Diameter 1.6 mm
Length 100 mm



Silverline®

Scientific Information

Problem

Infections are a serious complication of external ventricular drainage. Even with best surgical technique and utmost aseptic handling measures infections cannot be avoided completely (1, 3, 5).

For ventricular catheters impregnated with antibiotics a reduction of the rate of infection was described (7). In the product used for that study a combination of two antibiotics is used (Rifampicin and Minoocyclin). Also a product is on the market that uses a combination of Rifampicin and Clindamycin. However, antibiotics can generate resistancies when used prophylactically. Also an increasing proportion of device-related infections are being caused by *Candida* (9).

Method

Silverline Catheters are impregnated with 2% of nanoparticles of silver and an unsoluble silver salt. Silver ions are a strong antiseptic with a broad spectrum (4, 6). Should a resistancy to silver develop under the influence of the the silver ions, the antibiotic therapy will still be possible.

The nanometer-size particles have a large specific surface resulting in the release of silver ions in antiseptic concentration on the surface of the catheter.

The silver salt starts to release silver ions in antiseptic concentration immediately after placement of the catheter. The release of silver ions from the silver particles starts several hours later, when the ion release from the silver salt slows down. The release of silver ions from the silver particles then continues over a long interval.

Silver inactivates the enzyme systems with SH-groups of the cell wall that are active in the transmembraneous energy metabolism. Silver ions block the breathing chain of the bacteria (6).

Results

Silverline Ventricular Catheters were tested with the roll-culture method (2). Catheters were subjected to different bacteria suspensions. After that, roll-cultures were done in fixed time intervals. For control purposes non-antimicrobial catheters were subjected to the the same test.

Silverline Catheters could demonstrate their effect in a spectacular manner:

- Staph. aureus,
- Staph. epidermidis,
- Methicillin resistant Staphylococcus,
- E. coli,
- Pseudomonas aeruginosa,
- Candida albicans, and
- Candida glabrata

were eradicated (8).

The picture shows cultures of the Silverline Catheter and of a non-antimicrobial control made of the same material after subjecting them to a suspension of *Candida albicans* and *Candida glabrata*. The top rows are cultures of the Silverline Catheters after certain time intervals. The lower row is the control.

Discussion

The antimicrobial effect of Silverline Catheters could be shown in the in-vitro trial in a spectacular manner. The broad spectrum of efficacy against all hospital bacteria gives the product a very good effect.

Silver as an antimicrobial substance has been used in central venous catheters for some time (Erlangen Silver Catheter). A commercially available product that uses the same combination of antiseptic substances which is used in the Silverline Catheters, showed a significant reduction in the infection rate in controlled studies (2).

Because of its broad band effect with no danger of resistancies the Silverline Catheters are a new element for fighting infection of ventricular catheters.

Bibliography

1. Anneke, A.: Infektiöse Komplikationen von Liquoraußendrainagen. Dissertation Universität Heidelberg, 1999
2. Guggenbichler, J.P.; Juhl, G.; Braun, G.G.; Fraß, M.; Künstle, O.A.; Plötz, J.; Saffartzik, W.; Steinhäuser, M.; Wenisch, C.: Klinische Untersuchungen mit einem neuen Nano-Silber imprägnierten zentralvenösen Katheter. Hyg Med 28, 2003, 228-34
3. Lozier AP, Sciacca RR, Romagnoli MF, Connolly ES Ventriculostomy-related infections: A critical review of literature. Neurosurgery 51, 2002, 170-182
4. Naegeli, C. v.: Über die oligodynamischen Erscheinungen an lebenden Zellen. Neue Denkschr. Allg. Schweiz. Ges. Naturwiss. 33, 1893, 174 - 182
5. Pfisterer, W.; Mühlbauer, M.; Czech, T.; Reinprecht, A.: Early diagnosis of external ventricular drainage infection: results of a prospective study. J Neurol Neurosurg Psychiatry 74, 2003, 929-932
6. Thurmann, R.B.; Gerba, C.H.P.: The molecular mechanisms of copper and silver ion disinfection of bacteria and viruses. Crit. Rev. Environmental Control, 18, 1989, 295-315
7. Zabramski J.M.; Darouiche, R.O. et al.: Efficacy of antimicrobial-impregnated external ventricular catheters: a prospective, randomized, controlled trial. J Neurosurg, 98, 2003, 725-730
8. Zschaler, R.: Testing of the Antimicrobial Effect of Catheter Tubing with a Roll Culture Method. www.spiegelberg.de/home/documents/Zschaler.pdf
9. Kojic, E.M.; Darouiche, R.O.: Candida infections of medical devices. Clin. Microbiol. Rev., 17, 2004, 255-267

